

Phospho-BCL2 (S70) Antibody



PACO06930

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:**Gene ID:**

BCL2

Uniprot

P10415

Synonyms:

BCL2; Apoptosis regulator Bcl-2


Immunogen:

Synthesized peptide derived from human Bcl-2 around the phosphorylation site of S70.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Phospho-Bcl-2 (S70) Polyclonal Antibody.



Western Blot analysis of JK cells using Phospho-Bcl-2 (S70) Polyclonal Antibody.