

Phospho-PTK2 (S843) Antibody



PACO06974

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:**Gene ID:**

PTK2

Uniprot

Q05397

Synonyms:

PTK2; FAK; FAK1; Focal adhesion kinase 1; FADK 1; Focal adhesion kinase-related nonkinase; FRNK; Protein phosphatase 1 regulatory subunit 71; PPP1R71; Protein-tyrosine kinase 2; p125FAK; pp125FAK


Immunogen:

Synthesized peptide derived from human FAK around the phosphorylation site of S843.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of customer's lysis using Phospho-FAK (S843) antibody.



Western blot analysis of HELA KB 293T 3T3 lysis using Phospho-FAK (S843) antibody.