

USP44 Antibody



PACO06991

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

USP44

Uniprot

Q9H0E7

Synonyms:

USP44; Ubiquitin carboxyl-terminal hydrolase 44; Deubiquitinating enzyme 44; Ubiquitin thioesterase 44; Ubiquitin-specific-processing protease 44


Immunogen:

Synthesized peptide derived from the Internal region of human USP44.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of mouse-kidney lysis using USP44 antibody.



Western Blot analysis of 293T cells using USP44 Polyclonal Antibody.