

AGPAT9 Antibody



PACO06992

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

GPAT3

Uniprot

Q53EU6

Synonyms:

AGPAT9; GPAT3; MAG1; HMFN0839; Glycerol-3-phosphate acyltransferase 3; GPAT-3; 1-acylglycerol-3-phosphate O-acyltransferase 9; 1-AGP acyltransferase 9; 1-AGPAT 9; Acyl-CoA: glycerol-3-phosphate acyltransferase 3; hGPAT3; Lung cancer metastas


Immunogen:

Synthesized peptide derived from the C-terminal region of human LPAAT- theta.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of SH-SY5Y 3T3 lysis using LPAAT- theta antibody.



Western Blot analysis of HEPG2 using LPAAT- theta Polyclonal Antibody.