

# TRIM72 Antibody



PACO06998

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

WB, IHC

**Recommended dilutions:**

WB:1:1000, IHC:1:200-500

**Protein Background:**

**Gene ID:**

TRIM72

**Uniprot**

Q6ZMU5

**Synonyms:**

Tripartite motif-containing protein 72

**Immunogen:**


Recombinant Protein.

**Storage:**


PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of Mouse skeletal muscle using TRIM72 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Mouse Pancreas Tissue using TRIM72 Polyclonal Antibody.