

BEST2 Antibody



PACO07013

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

WB:1:500-1:1000, IHC:1:200-1:500

Protein Background:**Gene ID:**

BEST2

Uniprot

Q8NFU1

Synonyms:

Bestrophin-2


Immunogen:

Synthetic peptide.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue using BEST2 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using BEST2 Polyclonal Antibody.