

STAT1 Antibody



PACO07046

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

WB:1:1000-1:2000, IHC:1:200-1:500

Protein Background:**Gene ID:**

STAT1

Uniprot

P42224

Synonyms:

Signal transducer and activator of transcription 1-alpha/beta


Immunogen:

Recombinant Protein.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) HeLa, 2) 3T3, 3) Rat Heart Tissue using Stat1 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Stat1 Polyclonal Antibody.