

PDE4B/PDE4C/PDE4D Antibody



PACO07067

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

PDE4B/PDE4C/PDE4D

Uniprot

Q07343/Q08493/Q08499

Synonyms:

PDE4B; DPDE4; cAMP-specific 3', 5'-cyclic phosphodiesterase 4B; DPDE4; PDE32; PDE4C; DPDE1; cAMP-specific 3',5'-cyclic phosphodiesterase 4C; DPDE1; PDE21; PDE4D; DPDE3; cAMP-specific 3',5'-cyclic phosphodiesterase 4D; DPDE3; PDE43


Immunogen:

Synthesized peptide derived from human PDE4B/C/D around the non-phosphorylation site of S133/119/190.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using PDE4B/C/D Polyclonal Antibody.



Western Blot analysis of K562 cells using PDE4B/C/D Polyclonal Antibody.