

ACTB Antibody



PACO07078

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:1000-1:4000,
IHC:1:100-1:300

Protein Background:

Gene ID:

ACTB

Uniprot

P60709

Synonyms:

ACTB; Actin; cytoplasmic 1; beta-actin


Immunogen:

Synthesized peptide derived from the N-terminal region of human Actin β .


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Western Blot analysis of HELA cells using Actin beta Polyclonal Antibody.



Western blot analysis of customer's lysis using Actin beta antibody.



Western Blot analysis of HepG2 cells using Actin beta Polyclonal Antibody.