

AKT1/AKT2/AKT3 Antibody



PACO07079

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

AKT1/AKT2/AKT3

Uniprot

P31749/P31751/Q9Y243

Synonyms:

AKT1; PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-Akt; RAC-PK-alpha; AKT2; RAC-beta serine/threonine-protein kinase; Protein kinase Akt-2; Protein kinase B


Immunogen:

Synthesized peptide derived from human Akt around the non-phosphorylation site of S473.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Western Blot analysis of NIH-3T3 cells using Akt Polyclonal Antibody.



Western Blot analysis of KB A549 using Akt Polyclonal Antibody.



Western Blot analysis of HeLa cells using Akt Polyclonal Antibody.