

NRG3 Antibody



PACO07126

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

NRG3

Uniprot

P56975

Synonyms:

NRG3; Pro-neuregulin-3, membrane-bound isoform; Pro-NRG3


Immunogen:

Synthesized peptide derived from the Internal region of human Neuregulin-3.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Western Blot analysis of mouse brain K562 cells using Neuregulin-3 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded mouse-brain, antibody was diluted at 1:200.



Western Blot analysis of MOUSE-BRAIN cells using Neuregulin-3 Polyclonal Antibody.