

ABL1 Antibody



PACO07133

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

ABL1

Uniprot

P00519

Synonyms:

ABL1; ABL; JTK7; Tyrosine-protein kinase ABL1; Abelson murine leukemia viral oncogene homolog 1; Abelson tyrosine-protein kinase 1; Proto-oncogene c-Abl; p150


Immunogen:

Synthesized peptide derived from human Abl1 around the non-phosphorylation site of T735.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Western Blot analysis of cos-7 cells using Abl1 Polyclonal Antibody.



Western blot analysis of HELA KB SH-SY5Y 293T 3T3 lysis using Abl1 antibody.



Western Blot analysis of COS cells using Abl1 Polyclonal Antibody.