

TUBA1A/TUBA1B Antibody



PACO07151

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

TUBA1A/TUBA1B

Uniprot

Q71U36/P68363

Synonyms:

TUBA1A; TUBA3; Tubulin alpha-1A chain; alpha-tubulin 3; Tubulin B-alpha-1; Tubulin alpha-3 chain; TUBA1B; Tubulin alpha-1B chain; alpha-tubulin ubiquitous; Tubulin K-alpha-1; Tubulin alpha-ubiquitous chain


Immunogen:

Synthesized peptide derived from the C-terminal region of human Tubulin α .


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Western Blot analysis of HeLa A549 cells using Tubulin alpha Polyclonal Antibody.



Western Blot analysis of HELA HEPG2 using Tubulin alpha Polyclonal Antibody.



Western Blot analysis of HepG2 cells using Tubulin alpha Polyclonal Antibody.