

CACNA1I Antibody




PACO07169

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse, Rat	Gene ID: CACNA1I
Source: Rabbit	Uniprot Q9P0X4
Isotype: IgG	Synonyms:
Applications: ELISA, IHC	Immunogen: Synthetic peptide.
Recommended dilutions: ELISA:1:2000-1:10000, IHC-p:1:100-1:200	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using Cav3.3Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Cav3.3Rabbit pAb diluted at 1:200.