

# KCNK4 Antibody



PACO07182

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC-p:1:100-1:200

**Protein Background:**

**Gene ID:**

KCNK4

**Uniprot**

Q9NYG8

**Synonyms:**

KCNK4; TRAAK; Potassium channel subfamily K member 4; TWIK-related arachidonic acid, stimulated potassium channel protein; TRAAK; Two pore potassium channel KT4.1; Two pore K(+) channel KT4.1

**Immunogen:**

Synthetic peptide.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Immunohistochemical analysis of paraffin-embedded Mouse BrainTissue using KCNK4 (TRAAK) Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Rat BrainTissue using KCNK4 (TRAAK) Rabbit pAb diluted at 1:200.