

ADAMTS2 Antibody



PACO07203

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

ADAMTS2

Uniprot

O95450

Synonyms:

ADAMTS2; PCINP; PCPNI; A disintegrin and metalloproteinase with thrombospondin motifs 2; ADAM-TS 2; ADAM-TS2; ADAMTS-2; Procollagen I N-proteinase; PC I-NP; Procollagen I/II amino propeptide-processing enzyme; Procollagen N-endopeptidase; pNPI


Immunogen:

Synthesized peptide derived from the C-terminal region of human ADAMTS-2. at AA range: 1140-1211.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of HepG2 K562 HeLa Colo using ADAMTS-2 antibody. Secondary antibody was diluted at 1:20000.