

ADAMTS18 Antibody



PACO07236

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IHC-
p:1:100-1:300

Protein Background:

Gene ID:

ADAMTS18

Uniprot

Q8TE60

Synonyms:

ADAMTS18; ADAMTS21; A disintegrin and metalloproteinase with thrombospondin motifs 18; ADAM-TS 18; ADAM-TS18; ADAMTS-18


Immunogen:

Synthesized peptide derived from the C-terminal region of human ADAMTS-18. at AA range: 1030-1110.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded human-skin, antibody was diluted at 1:200.



Western blot analysis of 3T3 using ADAMTS-18 antibody. Secondary antibody was diluted at 1:20000.