

GLB1 Antibody



PACO07238

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IHC-
p:1:100-1:300

Protein Background:**Gene ID:**

GLB1

Uniprot

P16278

Synonyms:

GLB1; ELNR1; beta-galactosidase; acid, beta-galactosidase; Lactase; Elastin receptor 1


Immunogen:

Synthesized peptide derived from the Internal region of human Galactosidase β .


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Immunohistochemical analysis of paraffin-embedded human-breast-cancer, antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded human-colont-cancer, antibody was diluted at 1:200.



Western blot analysis of 3T3 using Galactosidase beta antibody. Secondary antibody was diluted at 1:20000.