

# Acetyl-EPAS1 (K385) Antibody



PACO07245

---

## Product Information

<b>Size:</b> 50ug	<b>Protein Background:</b>
<b>Reactivity:</b> Human, Mouse, Rat	<b>Gene ID:</b> EPAS1
<b>Source:</b> Rabbit	<b>Uniprot</b> Q99814
<b>Isotype:</b> IgG	<b>Synonyms:</b> endothelial PAS domain protein 1
<b>Applications:</b> ELISA, WB	<b>Immunogen:</b> Synthetic peptide of Acetyl-EPAS-1 (K385).
<b>Recommended dilutions:</b> ELISA:1:10000-1:20000, WB:1:500-1:2000	<b>Storage:</b> Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A