

TRAF1 Antibody




PACO07264

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse	Gene ID: TRAF1
Source: Rabbit	Uniprot Q13077
Isotype: IgG	Synonyms: TNF receptor-associated factor 1 (Epstein-Barr virus-induced protein 6)
Applications: ELISA, WB	Immunogen: Synthesized peptide derived from TRAF1 at AA range: 191-240.
Recommended dilutions: ELISA:1:10000-1:20000, WB:1:500-1:2000	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HEPG2 cells using TRAF1 Polyclonal Antibody diluted at 1:500. Secondary antibody was diluted at 1:20000.