

FAM3D Antibody



PACO07266

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000-1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

FAM3D

Uniprot

Q96BQ1

Synonyms:

Protein FAM3D


Immunogen:

Synthesized peptide derived from FAM3D at AA range: 121-170.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using FAM3D Polyclonal Antibody diluted at 1:500. Secondary antibody was diluted at 1:20000.