

HSP90AA1 Antibody



PACO07302

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC-p:1:50-1:200

Protein Background:**Gene ID:**

HSP90AA1

Uniprot

P07900

Synonyms:

Heat shock protein HSP 90-alpha (Heat shock 86 kDa) (HSP 86) (HSP86) (Renal carcinoma antigen NY-REN-38)


Immunogen:

Recombinant Protein of HSP90 α .


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using HSP90 alpha Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using HSP90alpha Rabbit pAb diluted at 1:200.