

NFE2L2 Antibody



PACO07306

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC-p:1:50-1:200

Protein Background:**Gene ID:**

NFE2L2

Uniprot

Q16236

Synonyms:

Nuclear factor erythroid 2-related factor 2 (NF-E2-related factor 2) (NFE2-related factor 2) (HEBP1) (Nuclear factor, erythroid derived 2, like 2)


Immunogen:

Recombinant Protein of Nrf2.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using Nrf2 Rabbit pAb diluted at 1:500.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using Nrf2 Rabbit pAb diluted at 1:500.