

PTK2 Antibody



PACO07308

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC-p:1:50-1:200

Protein Background:

Gene ID:

PTK2

Uniprot

Q05397

Synonyms:

Focal adhesion kinase 1 (FADK 1) (EC 2.7.10.2) (Focal adhesion kinase-related nonkinase) (FRNK) (Protein phosphatase 1 regulatory subunit 71) (PPP1R71) (Protein-tyrosine kinase 2) (p125FAK) (pp125FAK)


Immunogen:

Recombinant Protein of FAK.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using FAK Rabbit pAb diluted at 1:500.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using FAK Rabbit pAb diluted at 1:500.