

# Acetyl-HNF4A/HNF4G (Lys127/79) Antibody



PACO07331

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:10000-1:20000, WB:1:500-1:2000

**Protein Background:****Gene ID:**

HNF4A/HNF4G

**Uniprot**

P41235/Q14541

**Synonyms:**

Hepatocyte nuclear factor 4-alpha (HNF-4-alpha) (Nuclear receptor subfamily 2 group A member 1) (Transcription factor 14) Hepatocyte nuclear factor 4-gamma (HNF-4-gamma) (Nuclear receptor subfamily 2 group A member 2)(TCF-14) (Transcription factor HNF-4)

**Immunogen:**


Synthetic Acetyl peptide from human protein at AA range: 127(HNF-4 &alpha;)/79(HNF-4 &gamma;).

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of mouse-lung mouse-brain mouse-heart Hela mouse-liver lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000.