

FABP2 Antibody



PACO07338

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000-1:20000, WB:1:500-1:2000,
IHC-p:1:50-1:300

Protein Background:**Gene ID:**

FABP2

Uniprot

P12104

Synonyms:

Fatty acid, binding protein, intestinal (Fatty acid, binding protein 2) (Intestinal-type fatty acid, binding protein) (I-FABP)


Immunogen:

Synthetic peptide from human protein at AA range: 90-132.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.


Product Images



Immunohistochemical analysis of paraffin-embedded human-small-intestine, antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded human-small-intestine, antibody was diluted at 1:200.



Western blot analysis of mouse-brain mouse-spinal-cord lysate, antibody was diluted at 2000. Secondary antibody was diluted at 1:20000.