

IL17B Antibody



PACO07385

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000-1:20000, IHC-p:1:50-1:300

Protein Background:

Gene ID:

IL17B

Uniprot

Q9UHF5

Synonyms:

Interleukin-17B (IL-17B) (Cytokine Zcyto7) (Interleukin-20) (IL-20) (Neuronal interleukin-17-related factor)


Immunogen:

Synthetic peptide from human protein at AA range: 1-50.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded human-small-intestine, antibody was diluted at 1:200.