

HLA-DMB Antibody



PACO07410

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000-1:20000, IHC-p:1:50-1:300

Protein Background:

Gene ID:

HLA-DMB

Uniprot

P28068

Synonyms:

HLA class II histocompatibility antigen, DM beta chain (MHC class II antigen DMB)
(Really interesting new gene 7 protein)


Immunogen:

Synthetic peptide from human protein at AA range: 40-100.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.


Product Images



Immunohistochemical analysis of paraffin-embedded Human-tonsil, antibody was diluted at 1:100.



Immunohistochemical analysis of paraffin-embedded Human-colon, antibody was diluted at 1:100.



Immunohistochemical analysis of paraffin-embedded Human-tonsil, antibody was diluted at 1:100.