

EREG Antibody




PACO07419

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse, Rat	Gene ID: EREG
Source: Rabbit	Uniprot O14944
Isotype: IgG	Synonyms: Proepiregulin [Cleaved into: Epiregulin (EPR)]
Applications: ELISA, IHC	Immunogen: Synthetic peptide from human protein at AA range: 71-120.
Recommended dilutions: ELISA:1:10000-1:20000, IHC-p:1:50-1:300	Storage: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded Human-liver, antibody was diluted at 1:100.