

EYA1/EYA4 Antibody



PACO07423

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000-1:20000, IHC-p:1:50-1:300

Protein Background:

Gene ID:

EYA1/EYA4

Uniprot

Q99502/O95677

Synonyms:


Immunogen:

Synthetic peptide from human protein at AA range: 271-320.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.


Product Images



Immunohistochemical analysis of paraffin-embedded Human-kidney, antibody was diluted at 1:100.



Immunohistochemical analysis of paraffin-embedded Human-brain, antibody was diluted at 1:100.



Immunohistochemical analysis of paraffin-embedded Human-kidney, antibody was diluted at 1:100.