

TPM1 Antibody



PACO07444

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000-1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

TPM1

Uniprot

P09493

Synonyms:

Tropomyosin alpha-1 chain (alpha-tropomyosin) (Tropomyosin-1)


Immunogen:

Synthesized peptide derived from human Tropomyosin α . at AA range: 101-150.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 1, mouse-lung 2, mouse-brain 3, mouse-spleen cells using primary antibody diluted at 1:500(4°C overnight). Secondary antibody: Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25°C, 1 hour).