

E2F1 Antibody



PACO07454

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:10000

Protein Background:**Gene ID:**

E2F1

Uniprot

Q01094

Synonyms:

E2F1 RBBP3


Immunogen:

Synthetic peptide from AA range: 100-170.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western Blot analysis of MOUSE-BRAIN cells using Antibody diluted at 500.