

FBXO7 Antibody



PACO07458

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:10000

Protein Background:

Gene ID:

FBXO7

Uniprot

Q9Y3I1

Synonyms:

FBXO7 FBX7


Immunogen:

Synthetic peptide from AA range: 371-420.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western Blot analysis of mouse-kidney cells using Antibody diluted at 1000.