

ATG5 Antibody



PACO07470

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:50-1:200

Protein Background:**Gene ID:**

ATG5

Uniprot

Q9H1Y0

Synonyms:

Autophagy protein 5 (APG5-like) (Apoptosis-specific protein)


Immunogen:

Recombinant Protein.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded Human Ovarian Carcinoma Tissue using ATG5 Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Hepatocarcinoma Tissue using ATG5 Rabbit pAb diluted at 1:200.