

ATG13 Antibody



PACO07472

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:50-1:200

Protein Background:

Gene ID:

ATG13

Uniprot

O75143

Synonyms:

Autophagy-related protein 13


Immunogen:

Recombinant Protein.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using ATG13 Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skeletal Muscle Tissue using ATG13 Rabbit pAb diluted at 1:200.