

Histone H2A. Z Antibody



PACO07480

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000-1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

H2AFZ

Uniprot

P0C0S5

Synonyms:

Histone H2A. Z (H2A/z)


Immunogen:

Synthetic peptide from AA range: 1-30.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of HELA Cell Lysate using antibody.