

TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA8/TUBA4A Antibody

PACO07481



Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000-1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA4A/TUBA8

Uniprot

Q71U36/P68363/Q9BQE3/Q13748/P68366/Q9NY65

Synonyms:

alpha-tubulin 1 antibodyALS22 antibodyB alpha 1 antibodybA408E5.3 antibodyH2 alpha antibodyHum a tub1 antibodyHum a tub2 antibodyLIS3 antibodyMGC171407 antibodyMGC55332 antibodyTBA4A_HUMAN antibodyTestis-specific alpha-tubulin antibodyTUBA1 antibodyTUBA1A antibodytuba1l antibodyTuba4a antibodyTubulin alpha


Immunogen:

Synthetic peptide from AA range: 380-420.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of HELA 293T Cell Lysate using antibody.