

MACC1 Antibody



PACO07500

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Acts as a transcription activator for MET and as a key regulator of HGF-MET signaling. Promotes cell motility, proliferation and hepatocyte growth factor (HGF)-dependent scattering in vitro and tumor growth and metastasis in vivo.

Gene ID:

MACC1

Uniprot

Q6ZN28

Synonyms:

SH3 domain-containing protein 7a5, MACC1

Immunogen:

Recombinant Human Metastasis-associated in colon cancer protein 1 protein.

Storage:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Product Images

N/A

N/A