

PACO07506

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Reversible hydration of carbon dioxide. Participates in pH regulation. May be involved in the control of cell proliferation and transformation. Appears to be a novel specific biomarker for a cervical neoplasia.

Gene ID:

CA9

Uniprot

Q16790

Synonyms:

Carbonate dehydratase IX, Carbonic anhydrase IX, CA-IX, CAIX, Membrane antigen MN, Renal cell carcinoma-associated antigen G250, CA9

Immunogen:

recombinant Carbonic anhydrase 9 protein.

Storage:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Product Images

N/A

N/A