

# A1BG Antibody



PACO07507

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

A1BG

**Uniprot**

P04217

**Synonyms:**

alpha-1-B glycoprotein; A1BG; A1B; ABG; DKFZp686F0970; GAB; HYST2477 ;

**Immunogen:**

Human A1BG.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A