

PACO07509

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

AAAS

Uniprot

Q9NRG9

Synonyms:

achalasia, adrenocortical insufficiency, alacrimia (Allgrove, triple-A); AAAS; AAA; AAASb; ADRACALA; ADRACALIN; ALADIN; DKFZp586G1624; GL003 ;

Immunogen:

Human AAAS.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A