

AADAT Antibody



PACO07511

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

AADAT

Uniprot

Q8N5Z0

Synonyms:

aminoadipate aminotransferase; AADAT; KAT2; KATII ;

Immunogen:

Human AADAT.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A