

AARSD1 Antibody



PACO07514

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

AARSD1

Uniprot

Q9BTE6

Synonyms:

alanyl-tRNA synthetase domain containing 1; AARSD1; MGC2744 ;

Immunogen:

Human AARSD1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A