

AATF Antibody



PACO07515

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

AATF

Uniprot

Q9NY61

Synonyms:

apoptosis antagonizing transcription factor; AATF; CHE-1; CHE1; DED ; Traube

Immunogen:

Human AATF.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A