

# ABCB9 Antibody



PACO07521

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

ABCB9

**Uniprot**

Q9NP78

**Synonyms:**

ATP-binding cassette, sub-family B (MDR/TAP), member 9; ABCB9; EST122234; KIAA1520; TAPL ;

**Immunogen:**

Human ABCB9.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A