

ABCC2 Antibody



PACO07522

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ABCC2

Uniprot

Q92887

Synonyms:

ATP-binding cassette, sub-family C (CFTR/MRP), member 2; ABCC2; ABC30; CMOAT; DJS; KIAA1010; MRP2; cMRP ;

Immunogen:

Human ABCC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A