

# ABHD14B Antibody



PACO07530

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ABHD14B

**Uniprot**

Q96IU4

**Synonyms:**

abhydrolase domain containing 14B; ABHD14B; CIB; MGC15429 ;

**Immunogen:**

Human ABHD14B.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A