

ABI1 Antibody



PACO07533

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

Protein Background:

Gene ID:

ABI1

Uniprot

Q8IZP0

Synonyms:

507aa,55kd;

Immunogen:

Human ABI1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A