

ABI2 Antibody



PACO07534

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ABI2

Uniprot

Q9NYB9

Synonyms:

abl interactor 2; ABI2; ABI-2; ABI2B; AIP-1; AbIBP3; SSH3BP2; argBPIA; argBPIB ;

Immunogen:

Human ABI2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A