

ABL2 Antibody



PACO07535

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ABL2

Uniprot

P42684

Synonyms:

v-abl Abelson murine leukemia viral oncogene homolog 2 (arg, Abelson-related gene); ABL2; ABLL; ARG; FLJ22224; FLJ31718; FLJ41441 ;

Immunogen:

Human ABL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A